



Distribution Specialists, Inc.
43671 Trade Center Place, Suite 110
Sterling, VA 20166-2118
Phone: 703.996.8142
Toll free: 726.8853
Fax: 703.996.8143
E-mail: sales@DSI-Ltd.com

PRODUCT FACT SHEET

Accessories for estimating ISO Cleanliness levels of Hydraulic Fluids and Lubricating Oils

100-power Microscope with reticles

Pall® Comparator Guide

Order code DS K30



This set of accessories will enable you to estimate the types and amounts of solid contaminants in a hydraulic oil or fluid lubricant as defined by the ISO Cleanliness Codes.


The graduated reticles of the illuminated 100-power microscope display the test patch in the same scale as the reference patch photos the Comparator Guide. The examples in the guide were prepared using actual field samples and laboratory equipment similar to that in this kit. The particle count information shown in the Guide, for reference only, is based on the revised ISO 4406 (FDIS) procedure.

With experience, users can typically estimate the contamination level of a test sample within one ISO Cleanliness Code.

Typically used with:



DSI Hydraulic Test Kit
Order code DS K2



DSI Hydraulic Particles Test Kit
Order code DS K25

